

For Research Use Only

# Cytokeratin 1-specific Recombinant antibody, PBS Only (Capture)

Catalog Number: 82972-3-PBS



## Basic Information

Catalog Number:

82972-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_006121

GeneID (NCBI):

3848

UNIPROT ID:

P04264

Full Name:

keratin 1

Calculated MW:

66 kDa

Purification Method:

Protein A purification

CloneNo.:

230261E1

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sandwich ELISA standard curve of MP00134-1, Cytokeratin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82972-3-PBS. Detection antibody: 82972-2-PBS. Standard: SY00109. Range: 0.313-20 ng/mL.