

For Research Use Only

CTGF Recombinant antibody, PBS Only (Detector)

Catalog Number: 82976-2-PBS

1 Publications



Basic Information

Catalog Number:

82976-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001901

GeneID (NCBI):

1490

UNIPROT ID:

P29279

Full Name:

connective tissue growth factor

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

230324C3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|----------|-----------|-------------------------|-------------|
| Guang Xu | 38946862 | Anal Cell Pathol (Amst) | WB |

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve ofMP00091-1, CTGF
Recombinant Matched Antibody Pair, PBS Only.
Capture antibody: 82976-4-PBS. Detection
antibody: 82976-2-PBS. Standard: SY01606. Range:
0.313-10 ng/mL