

For Research Use Only

ARID1B Recombinant antibody, PBS Only (Capture)



Catalog Number: 82979-2-PBS

Basic Information

Catalog Number:

82979-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30643

GenBank Accession Number:

NM_017519

GeneID (NCBI):

57492

UNIPROT ID:

Q8NFD5

Full Name:

AT rich interactive domain 1B (SWI1-like)

Calculated MW:

236 kDa

Purification Method:

Protein A purification

CloneNo.:

230222H3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00103-2, ARID1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82979-2-PBS. Detection antibody: 82979-5-PBS. Standard: Ag30643. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00103-3, ARID1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82979-2-PBS. Detection antibody: 82979-3-PBS. Standard: Ag30643. Range: 0.625-20 ng/mL