

For Research Use Only

ARID1B Recombinant antibody, PBS Only (Detector)



Catalog Number: 82979-3-PBS

Basic Information

Catalog Number:	82979-3-PBS	GenBank Accession Number:	NM_017519	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	57492	CloneNo.:	23022287
Source:	Rabbit	UNIPROT ID:	Q8NFD5		
Isotype:	IgG	Full Name:	AT rich interactive domain 1B (SWI1-like)		
Immunogen Catalog Number:	AG30643	Calculated MW:	236 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00103-3, ARID1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82979-2-PBS. Detection antibody: 82979-3-PBS. Standard: Ag30643. Range: 0.625-20 ng/mL