

For Research Use Only

# POU4F3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82980-5-PBS

## Basic Information

Catalog Number: 82980-5-PBS	GenBank Accession Number: BC104923	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5459	CloneNo.: 230345F9
Source: Rabbit	UNIPROT ID: Q15319	
Isotype: IgG	Full Name: POU class 4 homeobox 3	
Immunogen Catalog Number: AG16004	Calculated MW: 338 aa, 37 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve ofMP00087-2, POU4F3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82980-5-PBS. Detection antibody: 82980-2-PBS. Standard: Ag16004. Range: 0.625-20 ng/mL.