

For Research Use Only

# OGDHL Recombinant antibody, PBS Only (Detector)



Catalog Number: 82981-2-PBS

## Basic Information

Catalog Number: 82981-2-PBS	GenBank Accession Number: BC026320	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 55753	CloneNo.: 230228F2
Source: Rabbit	UNIPROT ID: Q9ULDO	
Isotype: IgG	Full Name: oxoglutarate dehydrogenase-like	
Immunogen Catalog Number: AG10794	Calculated MW: 1010 aa, 115 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00023-1, OGDHL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82981-3-PBS. Detection antibody: 82981-2-PBS. Standard: Ag10794. Range: 1.25-20 ng/mL.