

For Research Use Only

OGDHL Recombinant antibody, PBS Only (Capture)



Catalog Number: 82981-5-PBS

Basic Information

| | | |
|--------------------------------------|---|--|
| Catalog Number: 82981-5-PBS | GenBank Accession Number: BC026320 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 55753 | CloneNo.: 230228H6 |
| Source: Rabbit | UNIPROT ID: Q9ULDO | |
| Isotype: IgG | Full Name: oxoglutarate dehydrogenase-like | |
| Immunogen Catalog Number: AG10794 | Calculated MW: 1010 aa, 115 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00023-2, OGDHL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82981-5-PBS. Detection antibody: 82981-4-PBS. Standard: Ag10794. Range: 2.5-20 ng/mL.