

For Research Use Only

AFP Recombinant antibody, PBS Only (Detector)

Catalog Number: 82982-2-PBS



Basic Information

Catalog Number:

82982-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC027881

GeneID (NCBI):

174

UNIPROT ID:

P02771

Full Name:

alpha-fetoprotein

Calculated MW:

69 kDa

Purification Method:

Protein A purification

CloneNo.:

230435F11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00001-1, AFP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82982-1-PBS. Detection antibody: 82982-2-PBS. Standard: Eg0539. Range: 1.25-20 ng/mL



Cytometric bead array standard curve of MP00001-2, AFP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82982-3-PBS. Detection antibody: 82982-2-PBS. Standard: Eg0539. Range: 1.25-20 ng/mL



Cytometric bead array standard curve of MP00001-3, AFP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82982-4-PBS. Detection antibody: 82982-2-PBS. Standard: Eg0539. Range: 0.625-20 ng/mL