

For Research Use Only

TMEM132E Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 82983-1-PBS

Basic Information

Catalog Number:

82983-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26772

GenBank Accession Number:

NM_001304438

GeneID (NCBI):

124842

Full Name:

transmembrane protein 132E

Calculated MW:

116 kDa

Purification Method:

Protein A purification

CloneNo.:

230525H8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00054-2, TMEM132E Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82983-1-PBS. Detection antibody: 82983-2-PBS. Standard: Ag26772. Range: 1.25-20 ng/mL



Cytometric bead array standard curve ofMP00054-3, TMEM132E Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82983-4-PBS. Detection antibody: 82983-1-PBS. Standard: Ag26772. Range: 0.625-20 ng/mL