

For Research Use Only

ADAMDEC1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82985-3-PBS



Basic Information

Catalog Number:

82985-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12265

GenBank Accession Number:

BC074877

GeneID (NCBI):

27299

UNIPROT ID:

O15204

Full Name:

ADAM-like, decysin 1

Calculated MW:

470 aa, 53 kDa

Purification Method:

Protein A purification

CloneNo.:

230367B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00007-2, ADAMDEC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82985-5-PBS. Detection antibody: 82985-3-PBS. Standard:Ag12265. Range: 0.625-20 ng/mL.