

For Research Use Only

# Mouse IL-17F Recombinant antibody, PBS Only (Detector)

Catalog Number: 82987-3-PBS



## Basic Information

Catalog Number:

82987-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_145856.2

GeneID (NCBI):

257630

UNIPROT ID:

Q7TNI7-1

Full Name:

interleukin 17F

Calculated MW:

17KD

Purification Method:

Protein A purification

CloneNo.:

230027A9

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sandwich ELISA standard curve of MP01107-1, Mouse IL-17F Recombinant Matched Antibody Pair - PBS only. 82987-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0671. 82987-3-PBS was HRP conjugated as the detection antibody. Range: 3.9-125 pg/mL.