

For Research Use Only

# SLC25A3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82988-6-PBS



## Basic Information

Catalog Number:

82988-6-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0667

GenBank Accession Number:

BC003504

GeneID (NCBI):

5250

UNIPROT ID:

Q00325

Full Name:

solute carrier family 25  
(mitochondrial carrier; phosphate  
carrier), member 3

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

230054H7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00046-3, SLC25A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82988-6-PBS. Detection antibody: 82988-5-PBS. Standard: Ag0667. Range: 0.625-20 ng/mL