

For Research Use Only

# Mouse Hif-1 alpha Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 82989-3-PBS



## Basic Information

Catalog Number:

82989-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33620

GenBank Accession Number:

NM\_001313919

GeneID (NCBI):

15251

UNIPROT ID:

Q61221

Full Name:

hypoxia inducible factor 1, alpha subunit

Calculated MW:

93 kDa

Purification Method:

Protein A purification

CloneNo.:

230101H5

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP00056-1, Mouse HIF-1 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82989-3-PBS. Detection antibody: 82989-2-PBS. Standard: Ag33620. Range: 0.156-10 ng/mL.



Sandwich ELISA standard curve of MP00056-2, Mouse HIF-1 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82989-1-PBS. Detection antibody: 82989-3-PBS. Standard: Ag33620. Range: 0.156-10 ng/mL.