

For Research Use Only

POLR2K Recombinant antibody, PBS Only (Detector)

Catalog Number: 82995-1-PBS



Basic Information

Catalog Number:

82995-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33765

GenBank Accession Number:

BC018157

GeneID (NCBI):

5440

UNIPROT ID:

P53803

Full Name:

polymerase (RNA) II (DNA directed)

polypeptide K, 7.0kDa

Purification Method:

Protein A purification

CloneNo.:

230147E8

Applications

Tested Applications:

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01314-1, POLR2K Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82995-3-PBS. Detection antibody: 82995-1-PBS. Standard: Ag33765. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP01314-2, POLR2K Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82995-2-PBS. Detection antibody: 82995-1-PBS. Standard: Ag33765. Range: 0.313-20 ng/mL