

For Research Use Only

FOXO1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82997-3-PBS



Basic Information

Catalog Number:

82997-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13296

GenBank Accession Number:

BC021981

GeneID (NCBI):

2308

UNIPROT ID:

Q12778

Full Name:

forkhead box O1

Calculated MW:

655 aa, 70 kDa

Purification Method:

Protein A purification

CloneNo.:

230183H6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00024-1, FOXO1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82997-1-PBS. Detection antibody: 82997-3-PBS. Standard: Ag13296. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00024-2, FOXO1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82997-2-PBS. Detection antibody: 82997-3-PBS. Standard: Ag13296. Range: 1.25-20 ng/mL