

For Research Use Only

SLC35E4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82999-1-PBS

Basic Information

| | | |
|--------------------------------------|---|--|
| Catalog Number: 82999-1-PBS | GenBank Accession Number: BC040191 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 339665 | CloneNo.: 240089A12 |
| Source: Rabbit | UNIPROT ID: Q6ICL7 | |
| Isotype: IgG | Full Name: solute carrier family 35, member E4 | |
| Immunogen Catalog Number: AG34360 | Calculated MW: 350 aa, 37 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00107-1, SLC35E4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82999-2-PBS. Detection antibody: 82999-1-PBS. Standard: Ag34360. Range: 1.25-20 ng/mL