

For Research Use Only

SLC35E4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82999-2-PBS

Basic Information

Catalog Number: 82999-2-PBS	GenBank Accession Number: BC040191	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 339665	CloneNo.: 240089C12
Source: Rabbit	UNIPROT ID: Q6ICL7	
Isotype: IgG	Full Name: solute carrier family 35, member E4	
Immunogen Catalog Number: AG34360	Calculated MW: 350 aa, 37 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00107-1, SLC35E4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82999-2-PBS. Detection antibody: 82999-1-PBS. Standard: Ag34360. Range: 1.25-20 ng/mL