

For Research Use Only

UNC13A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83000-1-PBS

Basic Information

Catalog Number:

83000-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001080421

GeneID (NCBI):

23025

UNIPROT ID:

Q9UPW8

Full Name:

unc-13 homolog A (C. elegans)

Calculated MW:

193 kDa

Purification Method:

Protein A purification

CloneNo.:

240017A4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00100-2, UNC13A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83000-4-PBS. Detection antibody: 83000-1-PBS. Standard: SY01073. Range: 0.156-5 ng/mL



Cytometric bead array standard curve ofMP00100-3, UNC13A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83000-2-PBS. Detection antibody: 83000-1-PBS. Standard: SY01073. Range: 0.156-2.5 ng/mL