

For Research Use Only

UNC13A Recombinant antibody, PBS Only



Catalog Number: 83000-6-PBS

Basic Information

Catalog Number:

83000-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001080421

GeneID (NCBI):

23025

UNIPROT ID:

Q9UPW8

Full Name:

unc-13 homolog A (C. elegans)

Calculated MW:

193 kDa

Purification Method:

Protein A purification

CloneNo.:

230469G2

Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed SH-SY5Y cells using UNC13A/Munc13-1 antibody (83000-6-RR, Clone: 230469G2) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 83000-6-PBS in a different storage buffer formulation.