

For Research Use Only

Mouse PCSK9 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:83004-1-PBS



Basic Information

Catalog Number: 83004-1-PBS	GenBank Accession Number: NM_153565.2	Purification Method: Protein A purification
Concentration: 1mg/ml	GeneID (NCBI): 100102	CloneNo.: 230536B3
Source: Rabbit	UNIPROT ID: Q80W65	
Isotype: IgG	Full Name: proprotein convertase subtilisin/kexin type 9	
Immunogen Catalog Number: EG0640	Calculated MW: 75KD	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00059-1, Mouse PCSK9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83004-1-PBS. Detection antibody: 83004-2-PBS. Standard: Eg0640. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00059-2, Mouse PCSK9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83004-2-PBS. Detection antibody: 83004-1-PBS. Standard: Eg0640. Range: 0.625-20 ng/mL.