

For Research Use Only

# IL-12/IL-23 p40 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83008-3-PBS



## Basic Information

Catalog Number:

83008-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC067498

GeneID (NCBI):

3593

UNIPROT ID:

P29460

Full Name:

interleukin 12B (natural killer cell  
stimulatory factor 2, cytotoxic  
lymphocyte maturation factor 2, p40)

Calculated MW:

328 aa, 37 kDa

Purification Method:

Protein A purification

CloneNo.:

240008G6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00035-1, IL12B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83008-3-PBS. Detection antibody: 83008-2-PBS. Standard: Eg0431. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00035-3, IL12B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83008-3-PBS. Detection antibody: 83008-1-PBS. Standard: Eg0431. Range: 0.625-20 ng/mL