

For Research Use Only

golgin 97 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83009-2-PBS

Basic Information

Catalog Number:	83009-2-PBS	GenBank Accession Number:	BC032853	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	2800	CloneNo.:	230523A6
Source:	Rabbit	UNIPROT ID:	Q92805		
Isotype:	IgG	Full Name:	golgi autoantigen, golgin subfamily a, 1		
Immunogen Catalog Number:	AG3335	Calculated MW:	767 aa, 88 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00031-2, GOLGA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83009-3-PBS. Detection antibody: 83009-2-PBS. Standard: Ag3335. Range: 1.25-20 ng/mL.