

For Research Use Only

COL6A3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83010-2-PBS



Basic Information

Catalog Number:

83010-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34752

GenBank Accession Number:

NM_004369

GeneID (NCBI):

1293

UNIPROT ID:

P12111

Full Name:

collagen, type VI, alpha 3

Calculated MW:

344 kDa

Purification Method:

Protein A purification

CloneNo.:

230540C10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00027-2, COL6A3-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83010-2-PBS. Detection antibody: 83010-4-PBS. Standard: Ag34752. Range: 0.625-20 ng/mL.