

For Research Use Only

GCSH Recombinant antibody, PBS Only



Catalog Number: 83012-5-PBS

Basic Information

Catalog Number:

83012-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10174

GenBank Accession Number:

BC000790

GeneID (NCBI):

2653

UNIPROT ID:

P23434

Full Name:

glycine cleavage system protein H
(aminomethyl carrier)

Calculated MW:

19 kDa

Purification Method:

Protein A purification

CloneNo.:

230505A3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 83012-5-RR (GCSH antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83012-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83012-5-RR against Human GCSH were performed. The affinity constant is below 1 pM.