

For Research Use Only

GPR87 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83013-1-PBS

Basic Information

| | | |
|--------------------------------------|---|--|
| Catalog Number: 83013-1-PBS | GenBank Accession Number: BC132888 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 53836 | CloneNo.: 230453E8 |
| Source: Rabbit | UNIPROT ID: Q9BY21 | |
| Isotype: IgG | Full Name: G protein-coupled receptor 87 | |
| Immunogen Catalog Number: AG23214 | | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00019-2, GPR87 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83013-1-PBS. Detection antibody: 83013-3-PBS. Standard: Ag23214. Range: 0.625-20 ng/mL