

For Research Use Only

GPR87 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83013-2-PBS

Basic Information

Catalog Number:

83013-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23214

GenBank Accession Number:

BC132888

GeneID (NCBI):

53836

UNIPROT ID:

Q9BY21

Full Name:

G protein-coupled receptor 87

Purification Method:

Protein A purification

CloneNo.:

230453H2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00019-1, GPR87 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83013-2-PBS. Detection antibody: 83013-3-PBS. Standard: Ag23214. Range: 0.625-20 ng/mL