

For Research Use Only

Mouse CD68 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83014-3-PBS

Basic Information

Catalog Number:

83014-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001291058.1

GeneID (NCBI):

12514

UNIPROT ID:

P31996-1

Full Name:

CD68 antigen

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

230504F10

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00015-1, Mouse CD68 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83014-3-PBS. Detection antibody: 83014-2-PBS. Standard: Eg0787. Range: 0.625-20 ng/mL.