

For Research Use Only

# Mouse CD68 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83014-4-PBS

## Basic Information

Catalog Number:

83014-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001291058.1

GeneID (NCBI):

12514

UNIPROT ID:

P31996-1

Full Name:

CD68 antigen

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

230504H3

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve ofMP00015-2, Mouse CD68 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83014-4-PBS. Detection antibody: 83014-2-PBS. Standard: Eg0787. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00015-3, Mouse CD68 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83014-4-PBS. Detection antibody: 83014-1-PBS. Standard: Eg0787. Range: 0.625-20 ng/mL