

For Research Use Only

FLT4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83015-4-PBS

Basic Information

Catalog Number:

83015-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002020

GeneID (NCBI):

2324

UNIPROT ID:

P35916

Full Name:

fms-related tyrosine kinase 4

Calculated MW:

153 kDa

Purification Method:

Protein A purification

CloneNo.:

240054C9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00009-1, FLT4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83015-2-PBS. Detection antibody: 83015-4-PBS. Standard: SY00928. Range: 0.313-20 ng/mL