

For Research Use Only

# Cadherin-19 Recombinant antibody

Catalog Number: 83019-1-RR



## Basic Information

Catalog Number:

83019-1-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2000

GenBank Accession Number:

BC015877

GeneID (NCBI):

28513

UNIPROT ID:

Q9H159

Full Name:

cadherin 19, type 2

Calculated MW:

55 kDa; 87 kDa

Purification Method:

Protein A purification

CloneNo.:

230271D2

## Applications

Tested Applications:

ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

