

For Research Use Only

Cadherin-19 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83019-3-PBS



Basic Information

Catalog Number:

83019-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2000

GenBank Accession Number:

BC015877

GeneID (NCBI):

28513

UNIPROT ID:

Q9H159

Full Name:

cadherin 19, type 2

Calculated MW:

55 kDa; 87 kDa

Purification Method:

Protein A purification

CloneNo.:

230271D7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00044-1, Cadherin-19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83019-4-PBS. Detection antibody: 83019-3-PBS. Standard: Ag2000. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00044-2, Cadherin-19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83019-5-PBS. Detection antibody: 83019-3-PBS. Standard: Ag2000. Range: 0.625-20 ng/mL