

For Research Use Only

FDFT1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83020-3-PBS



Basic Information

Catalog Number:

83020-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3743

GenBank Accession Number:

BC029641

GeneID (NCBI):

2222

UNIPROT ID:

P37268

Full Name:

farnesyl-diphosphate

farnesyltransferase 1

Calculated MW:

417 aa, 48 kDa

Purification Method:

Protein A purification

CloneNo.:

230362F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00097-3, SQ5/FDFT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83020-3-PBS. Detection antibody: 83020-2-PBS. Standard: Ag3743. Range: 0.625-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83020-3-PBS against Human FDFT1 was performed. The affinity constant is below 1 pM.