

For Research Use Only

ACVR1B Recombinant antibody, PBS Only (Capture)

Catalog Number: 83026-1-PBS



Basic Information

Catalog Number:

83026-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0129

GenBank Accession Number:

BC000254

GeneID (NCBI):

91

UNIPROT ID:

P36896

Full Name:

activin A receptor, type IB

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

230260C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00098-2, ACVR1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83026-1-PBS. Detection antibody: 83026-4-PBS. Standard: Ag0129. Range: 0.313-20 ng/mL