

For Research Use Only

CLDN18 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83028-1-PBS



Basic Information

Catalog Number:

83028-1-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34049

GenBank Accession Number:

NM_001002026

GeneID (NCBI):

51208

UNIPROT ID:

P56856-2

Full Name:

claudin 18

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

230198A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00105-1, Claudin 18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83028-2-PBS. Detection antibody: 83028-1-PBS. Standard: Ag34049. Range: 0.313-20 ng/mL.



Cytometric bead array standard curve of MP00105-2, Claudin 18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83028-1-PBS. Detection antibody: 83028-3-PBS. Standard: Ag34049. Range: 0.313-20 ng/mL.