

For Research Use Only

BMPR1B Recombinant antibody, PBS Only (Capture)

Catalog Number: 83029-7-PBS



Basic Information

Catalog Number: 83029-7-PBS	GenBank Accession Number: BC047773	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 658	CloneNo.: 230244H7
Source: Rabbit	UNIPROT ID: O00238	Affinity: $K_D = 9.83 \times 10^{-10}$
Isotype: IgG	Full Name: bone morphogenetic protein receptor, type IB	$K_{off} = 4.89 \times 10^{-4}$ $K_{on} = 4.97 \times 10^5$
Immunogen Catalog Number: AG24902	Calculated MW: 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00108-1, BMPR1B Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83029-7-PBS. Detection antibody: 83029-2-PBS. Standard: Ag24902. Range: 0.313-10 ng/mL