

For Research Use Only

GABRA4 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83030-2-PBS

Basic Information

Catalog Number: 83030-2-PBS	GenBank Accession Number: BC035055	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2557	CloneNo.: 230302E1
Source: Rabbit	UNIPROT ID: P48169	
Isotype: IgG	Full Name: gamma-aminobutyric acid (GABA) A receptor, alpha 4	
Immunogen Catalog Number: AG3600	Calculated MW: 554 aa, 62 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00115-1, GABRA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83030-2-PBS. Detection antibody: 83030-1-PBS. Standard: Ag3600. Range: 0.313-20 ng/mL



Cytometric bead array standard curve ofMP00115-2, GABRA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83030-3-PBS. Detection antibody: 83030-2-PBS. Standard: Ag3600. Range: 0.625-20 ng/mL