

For Research Use Only

# GABRA4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83030-3-PBS

## Basic Information

Catalog Number: 83030-3-PBS	GenBank Accession Number: BC035055	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2557	CloneNo.: 230302F5
Source: Rabbit	UNIPROT ID: P48169	
Isotype: IgG	Full Name: gamma-aminobutyric acid (GABA) A receptor, alpha 4	
Immunogen Catalog Number: AG3600	Calculated MW: 554 aa, 62 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00115-2, GABRA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83030-3-PBS. Detection antibody: 83030-2-PBS. Standard: Ag3600. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00115-3, GABRA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83030-3-PBS. Detection antibody: 83030-1-PBS. Standard: Ag3600. Range: 1.25-20 ng/mL