

For Research Use Only

COL15A1 Recombinant antibody

Catalog Number: 83031-3-RR



Basic Information

Catalog Number:

83031-3-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34644

GenBank Accession Number:

BC157094

GeneID (NCBI):

1306

UNIPROT ID:

P39059

Full Name:

collagen, type XV, alpha 1

Purification Method:

Protein A purification

CloneNo.:

230248F4

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83031-3-RR against Human COL15A1 were performed. The affinity constant is 0.819 nM.