

For Research Use Only

# COL15A1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83031-4-PBS



## Basic Information

Catalog Number:

83031-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34644

GenBank Accession Number:

BC157094

GeneID (NCBI):

1306

UNIPROT ID:

P39059

Full Name:

collagen, type XV, alpha 1

Purification Method:

Protein A purification

CloneNo.:

230248G2

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00133-1, COL15A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83031-4-PBS. Detection antibody: 83031-2-PBS. Standard: Ag34644. Range: 0.625-20 ng/mL



Sandwich ELISA standard curve of MP00133-2, COL15A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83031-4-PBS. Detection antibody: 83031-3-PBS. Standard: Ag34644. Range: 1.25-20 ng/mL