

For Research Use Only

FGF-21 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83033-2-PBS



Basic Information

Catalog Number:

83033-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC018404

GeneID (NCBI):

26291

UNIPROT ID:

Q9NSA1

Full Name:

fibroblast growth factor 21

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

230236C9

Affinity:

$K_D = 3.72 \times 10^{-9}$

$K_{off} = 1.49 \times 10^{-3}$

$K_{on} = 4.02 \times 10^5$

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C .

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00005-1, FGF21 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83033-1-PBS. Detection antibody: 83033-2-PBS. Standard: Eg0102. Range: 31.25-2000 pg/mL