

For Research Use Only

# FGF21 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83033-5-PBS



## Basic Information

Catalog Number:

83033-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0102

GenBank Accession Number:

BC018404

GeneID (NCBI):

26291

UNIPROT ID:

Q9NSA1

Full Name:

fibroblast growth factor 21

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

250540H7

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00005-2, Human FGF-21 Recombinant Matched Antibody Pair - PBS only. 83033-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ego102. 83033-4-PBS was HRP conjugated as the detection antibody. Range: 19.5-1250 pg/mL



Serum of ten individual healthy human donors was measured. The FGF-21 concentration of detected samples was determined to be 155.3 pg/mL with a range of 27.2-887.5 pg/mL