

For Research Use Only

# CD53 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83034-2-PBS

## Basic Information

Catalog Number: 83034-2-PBS	GenBank Accession Number: NM_000560	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 963	CloneNo.: 230188B7
Source: Rabbit	UNIPROT ID: P19397	
Isotype: IgG	Full Name: CD53 molecule	
Immunogen Catalog Number: AG28227	Calculated MW: 24 aa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00010-1, CD53 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83034-3-PBS. Detection antibody: 83034-2-PBS. Standard: Ag28227. Range: 0.313-20 ng/mL