

For Research Use Only

# EIF2C4 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83035-1-PBS

## Basic Information

Catalog Number: 83035-1-PBS	GenBank Accession Number: BC152450	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 192670	CloneNo.: 230265A5
Source: Rabbit	UNIPROT ID: Q9HCK5	
Isotype: IgG	Full Name: eukaryotic translation initiation factor 2C, 4	
Immunogen Catalog Number: AG19211	Calculated MW: 861 aa, 97 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00014-1, EIF2C4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83035-1-PBS. Detection antibody: 83035-2-PBS. Standard: Ag19211. Range: 0.313-20 ng/mL



Cytometric bead array standard curve ofMP00014-3, EIF2C4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83035-3-PBS. Detection antibody: 83035-1-PBS. Standard: Ag19211. Range: 0.313-20 ng/mL