

For Research Use Only

EIF2C4 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83035-3-PBS



Basic Information

Catalog Number:

83035-3-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19211

GenBank Accession Number:

BC152450

GeneID (NCBI):

192670

UNIPROT ID:

Q9HCK5

Full Name:

eukaryotic translation initiation factor 2C, 4

Calculated MW:

861 aa, 97 kDa

Purification Method:

Protein A purification

CloneNo.:

230265H8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00014-2, EIF2C4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83035-2-PBS. Detection antibody: 83035-3-PBS. Standard: Ag19211. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00014-3, EIF2C4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83035-3-PBS. Detection antibody: 83035-1-PBS. Standard: Ag19211. Range: 0.313-20 ng/mL