

For Research Use Only

# Acetyl-Histone H2A (Lys9)

## Recombinant antibody, PBS Only

### (Detector)

Catalog Number: **83041-3-PBS**

#### Basic Information

Catalog Number: 83041-3-PBS	GenBank Accession Number: BC093836	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3012	CloneNo.: 230349D4
Source: Rabbit	UNIPROT ID: P04908	
Isotype: IgG	Full Name: histone cluster 1, H2ae	
	Calculated MW: 14 kDa	

#### Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

#### Background Information

##### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00085-2, HIST1H2AE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83041-1-PBS. Detection antibody: 83041-3-PBS. Standard: SY02140Ac. Range: 0.156-5 ng/mL