

For Research Use Only

# NIACR1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83044-1-PBS



## Basic Information

Catalog Number:

83044-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12425

GenBank Accession Number:

BC027965

GeneID (NCBI):

338442

UNIPROT ID:

Q8TDS4

Full Name:

niacin receptor 1

Calculated MW:

363 aa, 42 kDa

Purification Method:

Protein A purification

CloneNo.:

230163B7

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00117-1, NIACR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83044-1-PBS. Detection antibody: 83044-4-PBS. Standard: Ag12425. Range: 0.625-20 ng/mL