

For Research Use Only

# NIACR1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83044-4-PBS

## Basic Information

Catalog Number: 83044-4-PBS	GenBank Accession Number: BC027965	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 338442	CloneNo.: 230163G1
Source: Rabbit	UNIPROT ID: Q8TDS4	
Isotype: IgG	Full Name: niacin receptor 1	
Immunogen Catalog Number: AG12425	Calculated MW: 363 aa, 42 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00117-1, NIACR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83044-1-PBS. Detection antibody: 83044-4-PBS. Standard: Ag12425. Range: 0.625-20 ng/mL