

For Research Use Only

TMEM20 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83045-1-PBS

Basic Information

Catalog Number: 83045-1-PBS	GenBank Accession Number: BC104814	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 159371	CloneNo.: 230298C11
Source: Rabbit	UNIPROT ID: Q2M3R5	
Isotype: IgG	Full Name: transmembrane protein 20	
Immunogen Catalog Number: AG34357	Calculated MW: 365 aa, 40 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00120-3, TMEM20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83045-3-PBS. Detection antibody: 83045-1-PBS. Standard: Ag34357. Range: 0.625-20 ng/mL.