

For Research Use Only

# TMEM20 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83045-4-PBS

## Basic Information

Catalog Number: 83045-4-PBS	GenBank Accession Number: BC104814	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 159371	CloneNo.: 230298F1
Source: Rabbit	UNIPROT ID: Q2M3R5	
Isotype: IgG	Full Name: transmembrane protein 20	
Immunogen Catalog Number: AG34357	Calculated MW: 365 aa, 40 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00120-1, TMEM20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83045-3-PBS. Detection antibody: 83045-4-PBS. Standard: Ag34357. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00120-2, TMEM20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83045-4-PBS. Detection antibody: 83045-2-PBS. Standard: Ag34357. Range: 0.625-20 ng/mL