

For Research Use Only

Claudin 5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83046-3-PBS



Basic Information

Catalog Number:

83046-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29044

GenBank Accession Number:

BC032363

GeneID (NCBI):

7122

UNIPROT ID:

O00501

Full Name:

claudin 5

Calculated MW:

218 aa, 23 kDa

Purification Method:

Protein A purification

CloneNo.:

230250E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00128-2, Claudin 5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83046-1-PBS. Detection antibody: 83046-3-PBS. Standard: Ag29044. Range: 0.313-20 ng/mL.